




Methyl tert-Butyl Ether CH3OC (CH3) 3

Part No.:166

Methyl tert-Butyl Ether CH3OC (CH3)3 No.166



Performance

Measuring range	10 to 200 ppm	200 to 660 ppm	
Number of pump strokes	2(200 mL)	1(100 mL)	
Correction factor	1	3.3	
Sampling time	3 min	1.5 min	

Detecting limit : 2 ppm (2 pump strokes)
Colour change : Yellow → Pale blue

Operating conditions : Temperature 0 to 35 °C (32 to 95 °F) correction used

Relative humidity 0 to 90 % correction not used

Relative standard deviation: 10 % (for 10 to 50 ppm), 5 % (for 50 to 200 ppm)

Tube quantity and number of tests per box: 10 tubes for 10 tests

Shelf life: 24 months

Reaction principle

 $CH_3OC(CH_3)_3 + Cr^{6+} + H_2SO_4 \rightarrow Cr^{3+}$

Possible coexisting substances and their interferences

Substance	Concentration	Interference	Changes colour by itself to
Methanol	≥ 3 times	+	Pale blue (≥ 5 ppm)
1-Propanol	≧ 1/1	+	Two layers
2-Propanol	≥ 1/1	+	Three layers
Aromatic hydrocarbons (Benzene)	≤ 30 ppm	No*	Pale brown (whole layer)
Aromatic hydrocarbons (Toluene)	≥ 1.0 ppm	+(Two layers, Unclear demarcation)	Blackish brown (≧ 0.5 ppm)
n-Hexane	≤ 10 ppm	No*	Pale brown (whole layer)

^{*}When Benzene coexists in a gas concentration of 30 ppm, or if n-Hexane coexists in a gas concentration of 10 ppm the reagent stains as two layers. The tube reading is unaffected.

Other substance measurable with this detector tube

Substance	Correction	No. of pump strokes	Measuring range		
Methyl isothiocyanate	Factor: 3.98	1	39.8 to 796 ppm		
	Factor : 8.83	1/2	88.3 to 1766 ppm		

Calibration gas generation

Diffusion tube method